

CHTOP antibody

Product Information

Catalog No.:	FNab01069
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


None

Immunogen information

Immunogen:	chromosome 1 open reading frame 77
Synonyms:	Chromatin target of PRMT1 protein Friend of PRMT1 protein Small arginine- and glycine-rich protein (SRAG) CHTOP C1orf77 FOP
Observed MW:	26-30 kDa
Uniprot ID :	Q9Y3Y2

Application

Reactivity:	Human, Mouse
Tested Application:	ELISA, WB
Recommended dilution:	WB: 1:500-1:2000
Image:	



mouse spleen tissue were subjected to SDS PAGE followed by western blot with FNab01069(C1orf77 antibody) at dilution of 1:500