

TCEANC2 antibody

Product Information

Catalog No.:	FNab01070
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


None

Immunogen information


Immunogen:	chromosome 1 open reading frame 83
Synonyms:	Transcription elongation factor A N-terminal and central domain-containing protein 2 TCEANC2 C1orf83
Observed MW:	24 kDa, 30 kDa
Uniprot ID :	Q96MN5

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IHC
Recommended dilution:	WB: 1:500-1:2000; IHC: 1:20-1:200
Image:	



Immunohistochemistry of paraffin-embedded human brain using FNab01070(C1orf83 antibody) at dilution of 1:50



human brain tissue were subjected to SDS PAGE followed by western blot with FNab01070(C1orf83 antibody) at dilution of 1:300