

LRRC71 antibody

Product Information

Catalog No.:	FNab01071
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background

None

Immunogen information

Immunogen:	chromosome 1 open reading frame 92
Synonyms:	Leucine-rich repeat-containing protein 71 LRRC71 C1orf92
Observed MW:	62 kDa
Uniprot ID :	Q8N4P6

Application

Reactivity:	Human
Tested Application:	ELISA, WB, IF
Recommended dilution:	WB: 1:200-1:1000; IF: 1:20-1:200
Image:	

HeLa cells were subjected to SDS PAGE followed by western blot with FNab01071(C1orf92 Antibody) at dilution of 1:300

