

C1QBP antibody

Product Information


Catalog No.:	FNab01073
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Immunogen information


Immunogen:	complement component 1, q subcomponent binding protein
Synonyms:	Complement component 1 Q subcomponent-binding protein, mitochondrial ASF/SF2-associated protein p32 Glycoprotein gC1qBP (C1qBP) Hyaluronan-binding protein 1 Mitochondrial matrix protein p32 gC1q-R protein p33 (SF2AP32) C1QBP GC1QBP HABP1 SF2P32
Observed MW:	32 kDa
Uniprot ID :	Q07021

Application

Reactivity:	Human, Mouse
Tested Application:	ELISA, WB, IHC, IF
Recommended dilution:	WB: 1:500-1:2000; IHC: 1:20-1:200; IF: 1:20-1:200
Image:	



Immunohistochemistry of paraffin-embedded human tonsillitis slide using FNab01073(C1QBP Antibody) at dilution of 1:50



K-562 cells were subjected to SDS PAGE followed by western blot with FNab01073(C1QBP Antibody) at dilution of 1:1000