

## DNAAF9 antibody

### Product Information

Catalog No.:	FNab01083
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

### Background


May act as an effector for ARL3.

### Immunogen information

Immunogen:	chromosome 20 open reading frame 194
Synonyms:	Dynein axonemal assembly factor 9 (DNAAF9) DNAAF9 C20orf194
Observed MW:	132 kDa
Uniprot ID :	Q5TEA3

### Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB
Recommended dilution: WB:	1:500-1:2000
Image:	



mouse testis tissue were subjected to SDS PAGE followed by western blot with FNab01083( C20orf194 Antibody) at dilution of 1:1000