

RTF2 antibody

Product Information

Catalog No.:	FNab01087
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


None

Immunogen information


Immunogen:	chromosome 20 open reading frame 43
Synonyms:	Replication termination factor 2 (RTF2) Replication termination factor 2 domain-containing protein 1 RTF2 C20orf43 RTFDC1
Observed MW:	40 kDa
Uniprot ID :	Q9BY42

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IP
Recommended dilution:	WB: 1:500-1:1000; IP: 1:200-1:1000
Image:	



IP Result of anti-C20orf43 (IP:FNab01087, 3ug;
Detection:FNab01087 1:500) with HeLa cells
lysate 3440ug.



MCF-7 cells were subjected to SDS PAGE
followed by western blot with
FNab01087(C20orf43 antibody) at dilution of
1:2000