

MIS18A antibody

Product Information

Catalog No.:	FNab01089
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


None

Immunogen information

Immunogen:	chromosome 21 open reading frame 45
Synonyms:	Protein Mis18-alpha FAPP1-associated protein 1 MIS18A C21orf45 C21orf46 FASP1
Observed MW:	105 kDa
Uniprot ID :	Q9NYP9

Application

Reactivity:	Human, Mouse
Tested Application:	ELISA, WB, IF
Recommended dilution:	WB: 1:200-1:1000; IF: 1:20-1:200
Image:	



HeLa cells were subjected to SDS PAGE followed by western blot with FNab01089(C21orf45 Antibody) at dilution of 1:300