

## XXYLT1 antibody

### Product Information

Catalog No.:	FNab01106
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

### Background


None

### Immunogen information


Immunogen:	chromosome 3 open reading frame 21
Synonyms:	Xyloside xylosyltransferase 1 UDP-xylose:alpha-xyloside alpha-1,3-xylosyltransferase XXYLT1 C3orf21
Observed MW:	44 kDa
Uniprot ID :	Q8NBI6

### Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IHC, IF
Recommended dilution:	WB: 1:500-1:2000; IHC: 1:20-1:200; IF: 1:20-1:200
Image:	



Immunohistochemistry of paraffin-embedded human prostate hyperplasia using FNab01106(C3orf21 antibody) at dilution of 1:50



HeLa cells were subjected to SDS PAGE followed by western blot with FNab01106(C3orf21 antibody) at dilution of 1:800