

C3orf38 antibody

Product Information

Catalog No.:	FNab01108
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


May be involved in apoptosis regulation.

Immunogen information

Immunogen:	chromosome 3 open reading frame 38
Synonyms:	Uncharacterized protein C3orf38 C3orf38
Observed MW:	35-37 kDa
Uniprot ID :	Q5JPI3

Application

Reactivity:	Human, Mouse
Tested Application:	ELISA, WB
Recommended dilution: WB:	1:500-1:2000
Image:	



L02 cells were subjected to SDS PAGE followed by western blot with FNab01108(C3orf38 Antibody) at dilution of 1:600