

POMGNT2 antibody

Product Information

Catalog No.:	FNab01109
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


None

Immunogen information


Immunogen:	chromosome 3 open reading frame 39
Synonyms:	Protein O-linked-mannose beta-1,4-N-acetylglucosaminyltransferase 2 (POMGnT2) Extracellular O-linked N-acetylglucosamine transferase-like Glycosyltransferase-like domain-containing protein 2 POMGNT2 AGO61 C3orf39 EOGTL GTDC2
Observed MW:	67-70 kDa
Uniprot ID :	Q8NAT1

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IP
Recommended dilution:	WB: 1:500-1:1000; IP: 1:200-1:1000
Image:	



IP Result of anti-C3orf39 (IP:FNab01109, 3ug;
Detection:FNab01109 1:400) with mouse skeletal
muscle tissue lysate 2500ug.



mouse skeletal muscle tissue were subjected to
SDS PAGE followed by western blot with
FNab01109(C3orf39 Antibody) at dilution of
1:600