

C4(gamma chain) antibody

Product Information


Catalog No.:	FNab01112
Size:	100µg
Form:	liquid
Purification:	Protein A+G purification
Purity:	≥95% as determined by SDS-PAGE
Host:	Mouse
Clonality:	monoclonal
Clone ID:	4A9
IsoType:	IgG1
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Immunogen information


Immunogen:	complement component 4A(Rodgers blood group)
Synonyms:	Complement C4-A Acidic complement C4 C3 and PZP-like alpha-2-macroglobulin domain-containing protein 2 Complement C4 beta chain Complement C4-A alpha chain C4a anaphylatoxin C4b-A C4d-A Complement C4 gamma chain C4A CO4 CPAMD2
Observed MW:	31 kDa
Uniprot ID :	P0C0L4

Application

Reactivity:	Human
Tested Application:	ELISA, WB, IHC
Recommended dilution:	WB: 1:2000-1:10000; IHC: 1:200-1:1000
Image:	



Immunohistochemistry of paraffin-embedded human liver tissue slide using FNab01112(C4 (gamma chain) antibody) at dilution of 1:800



human plasma tissue were subjected to SDS PAGE followed by western blot with FNab01112(C4 (gamma chain) antibody at dilution of 1:8000