

## MGARP antibody

### Product Information

Catalog No.:	FNab01117
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

### Background


None

### Immunogen information


Immunogen:	chromosome 4 open reading frame 49
Synonyms:	Protein MGARP Corneal endothelium-specific protein 1 (CESP-1) Hypoxia up-regulated mitochondrial movement regulator protein Mitochondria-localized glutamic acid-rich protein Ovary-specific acidic protein MGARP C4orf49 CESP1 HUMMR OSAP
Observed MW:	33-36 kDa
Uniprot ID :	Q8TDB4

### Application

Reactivity:	Human
Tested Application:	ELISA, IHC, WB
Recommended dilution:	WB: 1:500-1:2000; IHC: 1:20-1:200
Image:	



Immunohistochemistry of paraffin-embedded human testis slide using FNab01117(C4orf49 Antibody) at dilution of 1:50



SKOV-3 cells were subjected to SDS PAGE followed by western blot with FNab01117(C4orf49 antibody) at dilution of 1:1000