

## FAXDC2 antibody

### Product Information

Catalog No.:	FNab01121
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

### Background


None

### Immunogen information


Immunogen:	chromosome 5 open reading frame 4
Synonyms:	Fatty acid hydroxylase domain-containing protein 2 FAXDC2 C5orf4
Observed MW:	40 kDa
Uniprot ID :	Q96IV6

### Application

Reactivity:	Human
Tested Application:	ELISA, WB, IHC
Recommended dilution:	WB: 1:500-1:2000; IHC: 1:20-1:200
Image:	



Immunohistochemistry of paraffin-embedded human kidney tissue slide using FNab01121(C5orf4 Antibody) at dilution of 1:50



COLO 320 cells were subjected to SDS PAGE followed by western blot with FNab01121(C5orf4 Antibody) at dilution of 1:300