

CEP85L antibody

Product Information

Catalog No.:	FNab01129
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


None

Immunogen information


Immunogen:	chromosome 6 open reading frame 204
Synonyms:	Centrosomal protein of 85 kDa-like Serologically defined breast cancer antigen NY-BR-15 CEP85L C6orf204
Observed MW:	56 kDa
Uniprot ID :	Q5SZL2

Application

Reactivity:	Human
Tested Application:	ELISA, WB, IHC
Recommended dilution:	WB: 1:500-1:2000; IHC: 1:20-1:200
Image:	



Immunohistochemistry of paraffin-embedded human testis tissue slide using FNab01129(C6orf204 Antibody) at dilution of 1:50



PC-13 cells were subjected to SDS PAGE followed by western blot with FNab01129(C6orf204 Antibody) at dilution of 1:600