

ARMT1 antibody

Product Information

Catalog No.:	FNab01130
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


None

Immunogen information


Immunogen:	chromosome 6 open reading frame 211
Synonyms:	Damage-control phosphatase ARMT1 Acidic residue methyltransferase 1 Protein-glutamate O-methyltransferase Sugar phosphate phosphatase ARMT1 ARMT1 C6orf211
Observed MW:	50-51 kDa
Uniprot ID :	Q9H993

Application

Reactivity:	Human
Tested Application:	ELISA, WB, IHC
Recommended dilution:	WB: 1:500-1:2000; IHC: 1:20-1:200
Image:	



Immunohistochemistry of paraffin-embedded human testis tissue slide using FNab01130(C6orf211 Antibody) at dilution of 1:50



MCF-7 cells were subjected to SDS PAGE followed by western blot with FNab01130(C6orf211 Antibody) at dilution of 1:600