

## ERMARD antibody

### Product Information

Catalog No.:	FNab01131
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

### Background


None

### Immunogen information

Immunogen:	chromosome 6 open reading frame 70
Synonyms:	Endoplasmic reticulum membrane-associated RNA degradation protein (ER membrane-associated RNA degradation protein) ERMARD C6orf70
Observed MW:	78 kDa
Uniprot ID :	Q5T6L9

### Application

Reactivity:	Human
Tested Application:	ELISA, IHC
Recommended dilution:	IHC: 1:20-1:200
Image:	



Immunohistochemistry of paraffin-embedded human colon cancer tissue slide using FNab01131(C6orf70 Antibody) at dilution of 1:50