

## C7 antibody

### Product Information

Catalog No.:	FNab01132
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

### Background


Constituent of the membrane attack complex(MAC) that plays a key role in the innate and adaptive immune response by forming pores in the plasma membrane of target cells. C7 serves as a membrane anchor.

### Immunogen information

Immunogen:	complement component 7
Synonyms:	Complement component C7 C7
Observed MW:	105 kDa
Uniprot ID :	P10643

### Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB
Recommended dilution: WB:	1:500-1:2000
Image:	



mouse thymus tissue were subjected to SDS PAGE followed by western blot with FNab01132(C7 antibody) at dilution of 1:500