

CYREN antibody

Product Information

Catalog No.:	FNab01136
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


None

Immunogen information

Immunogen:	chromosome 7 open reading frame 49
Synonyms:	Cell cycle regulator of non-homologous end joining (Cell cycle regulator of NHEJ) Modulator of retrovirus infection homolog CYREN C7orf49 MRI
Observed MW:	17 kDa
Uniprot ID :	Q9BWK5

Application

Reactivity:	Human
Tested Application:	ELISA, IHC
Recommended dilution:	IHC: 1:10-1:100
Image:	



Immunohistochemistry of paraffin-embedded human kidney tissue slide using FNab01136(C7orf49 Antibody) at dilution of 1:50