

CPED1 antibody

Product Information

Catalog No.:	FNab01138
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


None

Immunogen information


Immunogen:	chromosome 7 open reading frame 58
Synonyms:	Cadherin-like and PC-esterase domain-containing protein 1 CPED1 C7orf58
Observed MW:	89 kDa
Uniprot ID :	A4D0V7

Application

Reactivity:	Human
Tested Application:	ELISA, WB, IHC, IF
Recommended dilution:	WB: 1:500-1:3000; IHC: 1:20-1:200; IF: 1:20-1:200
Image:	



Immunohistochemistry of paraffin-embedded human stomach using FNab01138(C7orf58 antibody) at dilution of 1:100



A375 cells were subjected to SDS PAGE followed by western blot with FNab01138(C7orf58 antibody) at dilution of 1:500