

## C9orf72 antibody

### Product Information

Catalog No.:	FNab01148
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

### Background


May play a role in endosomal trafficking and autophagy.

### Immunogen information


Immunogen:	chromosome 9 open reading frame 72
Synonyms:	Guanine nucleotide exchange factor C9orf72 C9orf72 DENND9 DENNL72
Observed MW:	25-30 kDa
Uniprot ID :	Q96LT7

### Application

Reactivity:	Human
Tested Application:	ELISA, IHC, WB
Recommended dilution:	WB: 1:500-1:2000; IHC: 1:20-1:200
Image:	



Immunohistochemistry of paraffin-embedded human gliomas tissue slide using FNab01148(C9orf72 Antibody) at dilution of 1:50



SH-SY5Y cells were subjected to SDS PAGE followed by western blot with FNab01148(C9orf72 Antibody) at dilution of 1:300