

CA9 antibody

Product Information

Catalog No.:	FNab01154
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


None

Immunogen information


Immunogen:	carbonic anhydrase IX
Synonyms:	Carbonic anhydrase 9 Carbonate dehydratase IX Carbonic anhydrase IX (CA-IX, CAIX) Membrane antigen MN P54/58N Renal cell carcinoma-associated antigen G250 (RCC-associated antigen G250) pMW1 CA9 G250 MN
Observed MW:	50 kDa
Uniprot ID :	Q16790

Application

Reactivity:	Human, Mouse
Tested Application:	ELISA, WB, IHC
Recommended dilution:	WB: 1:500-1:2000; IHC: 1:20-1:200
Image:	



Immunohistochemistry of paraffin-embedded human renal cell carcinoma tissue slide using FNab01154(CA9 Antibody) at dilution of 1:200 heat mediated antigen retrieved with Tris-EDTA buffer(pH9).



mouse brain tissue were subjected to SDS PAGE followed by western blot with FNab01154(CA IX antibody) at dilution of 1:300