

CA10 antibody

Product Information

Catalog No.:	FNab01156
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


Does not have a catalytic activity.

Immunogen information


Immunogen:	carbonic anhydrase X
Synonyms:	Carbonic anhydrase-related protein 10 Carbonic anhydrase-related protein X (CA-RP X, CARP X) Cerebral protein 15 CA10
Observed MW:	40 kDa
Uniprot ID :	Q9NS85

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, IHC, WB
Recommended dilution:	WB: 1:500-1:2000; IHC: 1:20-1:200
Image:	



Immunohistochemistry of paraffin-embedded human gliomas using FNab01156(CA10 antibody) at dilution of 1:100



Transfected HEK-293 cells were subjected to SDS PAGE followed by western blot with FNab01156(CA10 Antibody) at dilution of 1:700