

CA12 antibody

Product Information

Catalog No.:	FNab01158
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


Reversible hydration of carbon dioxide.

Immunogen information


Immunogen:	carbonic anhydrase XII
Synonyms:	Carbonic anhydrase 12 Carbonate dehydratase XII Carbonic anhydrase XII (CA-XII) Tumor antigen HOM-RCC-3.1.3 CA12
Observed MW:	44 kDa, 75-80 kDa
Uniprot ID :	O43570

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IHC, FC
Recommended dilution:	WB: 1:500-1:2000; IHC: 1:20-1:200
Image:	



Immunohistochemistry of paraffin-embedded human kidney using FNab01158(CA12 antibody) at dilution of 1:50



BxPC-3 cells were subjected to SDS PAGE followed by western blot with FNab01158(CA12 antibody) at dilution of 1:300