

CA5B antibody

Product Information

Catalog No.:	FNab01164
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


Reversible hydration of carbon dioxide.

Immunogen information

Immunogen:	carbonic anhydrase VB, mitochondrial
Synonyms:	Carbonic anhydrase 5B, mitochondrial Carbonate dehydratase VB Carbonic anhydrase VB (CA-VB) CA5B
Observed MW:	32 kDa
Uniprot ID :	Q9Y2D0

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB
Recommended dilution:	WB: 1:500-1:2000
Image:	



mouse testis tissue were subjected to SDS PAGE followed by western blot with FNab01164(CA5B antibody) at dilution of 1:4000