

COQ8A antibody

Product Information

Catalog No.:	FNab01169
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


None

Immunogen information

Immunogen:	chaperone, ABC1 activity of bc1 complex homolog(S. pombe)
Synonyms:	Atypical kinase COQ8A, mitochondrial Chaperone activity of bc1 complex-like (Chaperone-ABC1-like) Coenzyme Q protein 8A aarF domain-containing protein kinase 3 COQ8A ADCK3 CABC1
Observed MW:	60 kDa
Uniprot ID :	Q8NI60

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IF
Recommended dilution:	WB: 1:500-1:2000; IF: 1:20-1:200
Image:	



Hela cells were subjected to SDS PAGE followed by western blot with FNab01169(CABC1 antibody) at dilution of 1:2000