

## CADM3 antibody

### Product Information

Catalog No.:	FNab01191
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

### Background

Involved in the cell-cell adhesion. Has both calcium-independent homophilic cell-cell adhesion activity and calcium-independent heterophilic cell-cell adhesion activity with IGSF4, NECTIN1 and NECTIN3. Interaction with EPB41L1 may regulate structure or function of cell-cell junctions(By similarity).


### Immunogen information

Immunogen:	cell adhesion molecule 3
Synonyms:	Cell adhesion molecule 3 Brain immunoglobulin receptor Immunoglobulin superfamily member 4B (IGSF4B) Nectin-like protein 1 (NECL-1) Synaptic cell adhesion molecule 3 (SynCAM3) TSLC1-like protein 1 (TSL1) CADM3 IGSF4B NECL1 SYNCAM3 TSLL1
Observed MW:	38-50 kDa
Uniprot ID :	Q8N126

### Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB
Recommended dilution: WB:	1:500-1:2000

Image:



Jurkat cells were subjected to SDS PAGE followed by western blot with FNa01191(CADM3 antibody) at dilution of 1:800