

## CALD1 antibody

### Product Information

Catalog No.:	FNab01199
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

### Background


None

### Immunogen information


Immunogen:	caldesmon 1
Synonyms:	Caldesmon (CDM) CALD1 CAD CDM
Observed MW:	61-65 kDa, 70-80 kDa
Uniprot ID :	Q05682

### Application

Reactivity:	Human, Mouse
Tested Application:	ELISA, WB, IF, IHC
Recommended dilution:	WB: 1:500-1:2000; IHC: 1:20-1:200; IF: 1:20-1:200
Image:	



Immunohistochemistry of paraffin-embedded human colon tissue slide using FNab01199(CALD1 Antibody) at dilution of 1:50



mouse small intestine tissue were subjected to SDS PAGE followed by western blot with FNab01199(CALD1 antibody) at dilution of 1:500