

CAPN11 antibody

Product Information

Catalog No.:	FNab01212
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background

None

Immunogen information

Immunogen:	calpain 11
Synonyms:	Calpain-11 Calcium-activated neutral proteinase 11 (CANP 11) CAPN11
Observed MW:	81 kDa
Uniprot ID :	Q9UMQ6

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB
Recommended dilution: WB:	1:500-1:2000
Image:	

PC-3 cells were subjected to SDS PAGE followed by western blot with FNab01212(Calpain 11 antibody) at dilution of 1:1000

