

CAPN2 antibody

Product Information

Catalog No.:	FNab01213
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


None

Immunogen information


Immunogen:	calpain 2,(m/II) large subunit
Synonyms:	Calpain-2 catalytic subunit Calcium-activated neutral proteinase 2 (CANP 2) Calpain M-type Calpain large polypeptide L2 Calpain-2 large subunit Millimolar-calpain (M-calpain) CAPN2 CANPL2
Observed MW:	80 kDa
Uniprot ID :	P17655

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IHC, IF
Recommended dilution:	WB: 1:500-1:2000; IHC: 1:20-1:200; IF: 1:20-1:200
Image:	



Immunohistochemistry of paraffin-embedded human colon cancer using FNab01213(CAPN2 antibody) at dilution of 1:50



Y79 cells were subjected to SDS PAGE followed by western blot with FNab01213(CAPN2 antibody) at dilution of 1:200