

CAPN3 antibody

Product Information

Catalog No.:	FNab01214
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background

None

Immunogen information

Immunogen:	calpain 3,(p94)
Synonyms:	Calpain-3 Calcium-activated neutral proteinase 3 (CANP 3) Calpain L3 Calpain p94 Muscle-specific calcium-activated neutral protease 3 New calpain 1 (nCL-1) CAPN3 CANP3 CANPL3 NCL1
Observed MW:	100 kDa
Uniprot ID :	P20807

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA
Recommended dilution:	None