

## CAPN6 antibody

### Product Information

Catalog No.:	FNab01216
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

### Background


None

### Immunogen information


Immunogen:	calpain 6
Synonyms:	Calpain-6 Calpain-like protease X-linked Calpamodulin (CalpM) CAPN6 CALPM CANPX
Observed MW:	66-75 kDa
Uniprot ID :	Q9Y6Q1

### Application


Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, IHC, WB, IP
Recommended dilution:	WB: 1:500-1:2000; IP: 1:500-1:1000; IHC: 1:50-1:200
Image:	



Immunohistochemistry of paraffin-embedded human placenta using FNab01216(CAPN6 antibody) at dilution of 1:50



IP Result of anti-Calpain 6 (IP:FNab01216, 4ug; Detection:FNab01216 1:300) with human placenta tissue lysate 1520ug.



human placenta tissue were subjected to SDS PAGE followed by western blot with FNab01216(CAPN6 antibody) at dilution of 1:300