

CAPN9 antibody

Product Information

Catalog No.:	FNab01218
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


None

Immunogen information


Immunogen:	calpain 9
Synonyms:	Calpain-9 Digestive tract-specific calpain New calpain 4 (nCL-4) Protein CG36 CAPN9 NCL4
Observed MW:	70 kDa
Uniprot ID :	O14815

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IHC, IF
Recommended dilution:	WB: 1:500-1:2000; IHC: 1:20-1:200; IF: 1:10-1:100
Image:	



Immunohistochemistry of paraffin-embedded human stomach cancer using FNab01218(CAPN9 antibody) at dilution of 1:100



HT-29 cells were subjected to SDS PAGE followed by western blot with FNab01218(CAPN9 antibody) at dilution of 1:1000