

## CAPNS2 antibody

### Product Information

Catalog No.:	FNab01219
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

### Background


None

### Immunogen information

Immunogen:	calpain, small subunit 2
Synonyms:	Calpain small subunit 2 (CSS2) Calcium-dependent protease small subunit 2 CAPNS2
Observed MW:	28 kDa
Uniprot ID :	Q96L46

### Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB
Recommended dilution:	WB: 1:200-1:1000
Image:	



mouse skin tissue were subjected to SDS PAGE followed by western blot with FNab01219(Calpain S2 antibody) at dilution of 1:300