

CAMK2A-Specific antibody

Product Information

Catalog No.:	FNab01229
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


None

Immunogen information


Immunogen:	calcium/calmodulin-dependent protein kinase II alpha
Synonyms:	Calcium/calmodulin-dependent protein kinase type II subunit alpha (CaM kinase II subunit alpha, CaMK-II subunit alpha) CAMK2A CAMKA KIAA0968
Observed MW:	50 kDa
Uniprot ID :	Q9UQM7

Application

Reactivity:	Human
Tested Application:	ELISA, WB, IHC, IF
Recommended dilution:	WB: 1:500-1:2000; IHC: 1:20-1:200; IF: 1:20-1:200
Image:	



Immunohistochemistry of paraffin-embedded human brain tissue slide using FNab01229(CAMK2A-Specific Antibody) at dilution of 1:50



HEK-293 cells were subjected to SDS PAGE followed by western blot with FNab01229(CAMK2A-Specific antibody) at dilution of 1:1000