

CAMK2G antibody

Product Information

Catalog No.:	FNab01234
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background

None

Immunogen information


Immunogen:	calcium/calmodulin-dependent protein kinase II gamma
Synonyms:	Calcium/calmodulin-dependent protein kinase type II subunit gamma (CaM kinase II subunit gamma, CaMK-II subunit gamma) CAMK2G CAMK CAMK-II CAMKG
Observed MW:	59 kDa
Uniprot ID :	Q13555

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, IHC, IF, WB
Recommended dilution:	WB: 1:500-1:2000; IHC: 1:20-1:200; IF: 1:20-1:200
Image:	



Immunohistochemistry of paraffin-embedded human pancreas cancer using FNab01234(CAMK2G antibody) at dilution of 1:100



SMMC-7721 cells were subjected to SDS PAGE followed by western blot with FNab01234(CAMK2G Antibody) at dilution of 1:2000