

CAMKV antibody

Product Information

Catalog No.:	FNab01239
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


Does not appear to have detectable kinase activity.

Immunogen information

Immunogen:	CaM kinase-like vesicle-associated
Synonyms:	CaM kinase-like vesicle-associated protein CAMKV
Observed MW:	54-75 kDa
Uniprot ID :	Q8NCB2

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB
Recommended dilution: WB:	1:500-1:2000
Image:	



mouse brain tissue were subjected to SDS PAGE followed by western blot with FNab01239(CAMKV antibody) at dilution of 1:1200