

## CAP2 antibody

### Product Information

Catalog No.:	FNab01247
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

### Background


May have a regulatory bifunctional role.

### Immunogen information


Immunogen:	CAP, adenylate cyclase-associated protein, 2(yeast)
Synonyms:	Adenylyl cyclase-associated protein 2 (CAP 2) CAP2
Observed MW:	53 kDa
Uniprot ID :	P40123

### Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IHC
Recommended dilution:	WB: 1:500-1:2000; IHC: 1:20-1:200
Image:	



Immunohistochemistry of paraffin-embedded human liver cancer tissue slide using FNab01247(CAP2 Antibody) at dilution of 1:50



human kidney tissue were subjected to SDS PAGE followed by western blot with FNab01247(CAP2 antibody) at dilution of 1:400