

CAPG-Specific antibody

Product Information

Catalog No.:	FNab01249
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


None

Immunogen information


Immunogen:	capping protein(actin filament), gelsolin-like
Synonyms:	Macrophage-capping protein Actin regulatory protein CAP-G CAPG AFCP MCP
Observed MW:	43 kDa
Uniprot ID :	P40121

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IHC
Recommended dilution:	WB: 1:500-1:2000; IHC: 1:20-1:200
Image:	



Immunohistochemistry of paraffin-embedded human spleen using FNab01249(CAPG-Specific antibody) at dilution of 1:50



mouse lung tissue were subjected to SDS PAGE followed by western blot with FNab01249(CAPG-Specific antibody) at dilution of 1:300