

CPE antibody

Product Information

Catalog No.:	FNab01266
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


None

Immunogen information


Immunogen:	carboxypeptidase E
Synonyms:	Carboxypeptidase E (CPE) Carboxypeptidase H (CPH) Enkephalin convertase Prohormone-processing carboxypeptidase CPE CBPE
Observed MW:	50 kDa, 52 kDa
Uniprot ID :	P16870

Application


Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IHC, IP
Recommended dilution:	WB: 1:500-1:2000; IP: 1:500-1:1000; IHC: 1:50-1:200
Image:	



Immunohistochemistry of paraffin-embedded human pituitary using FNab01266(CPE antibody) at dilution of 1:50



IP Result of anti-Carboxypeptidase E (IP:FNab01266, 3ug; Detection:FNab01266 1:500) with mouse brain tissue lysate 3500ug.



mouse brain tissue were subjected to SDS PAGE followed by western blot with FNab01266(CPE antibody) at dilution of 1:600