

CARM1 antibody

Product Information


Catalog No.:	FNab01275
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Immunogen information


Immunogen:	coactivator-associated arginine methyltransferase 1
Synonyms:	Histone-arginine methyltransferase CARM1 Coactivator-associated arginine methyltransferase 1 Protein arginine N-methyltransferase 4 CARM1 PRMT4
Observed MW:	63-66 kDa, 55 kDa
Uniprot ID :	Q86X55

Application


Reactivity:	Human, Mouse, Rat, Monkey
Tested Application:	ELISA, IHC, WB, IF
Recommended dilution:	WB: 1:1000-1:5000; IHC: 1:50-1:500; IF: 1:20-1:200
Image:	



Immunohistochemistry of paraffin-embedded human breast cancer using FNab01275(CARM1 antibody) at dilution of 1:50



Immunofluorescent analysis of NIH/3T3 cells using FNab01275 (CARM1 antibody) at dilution of 1:50 and Alexa Fluor 488-conjugated Goat Anti-Rabbit IgG(H+L)



NIH/3T cells were subjected to SDS PAGE followed by western blot with FNab01275(CARM1 antibody) at dilution of 1:500