

CSNK1D antibody

Product Information

Catalog No.:	FNab01279
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


None

Immunogen information


Immunogen:	casein kinase 1, delta
Synonyms:	Casein kinase I isoform delta (CKI-delta, CKId) Tau-protein kinase CSNK1D CSNK1D HCKID
Observed MW:	47 kDa
Uniprot ID :	P48730

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IF, IP
Recommended dilution:	WB: 1:500-1:3000; IF: 1:20-1:200
Image:	



IP Result of anti-Casein Kinase 1 delta
(IP:FNab01279, 3ug; Detection:FNab01279
1:1000) with HeLa cells lysate 2000ug.



HeLa cells were subjected to SDS PAGE followed
by western blot with FNab01279(Casein Kinase 1
delta antibody) at dilution of 1:600