

## CASP2 antibody

### Product Information

Catalog No.:	FNab01288
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

### Background


None

### Immunogen information


Immunogen:	caspase 2, apoptosis-related cysteine peptidase
Synonyms:	Caspase-2 (CASP-2) Neural precursor cell expressed developmentally down-regulated protein 2 (NEDD-2) Protease ICH-1 Caspase-2 subunit p18 Caspase-2 subunit p13 Caspase-2 subunit p12 CASP2 ICH1 NEDD2
Observed MW:	50 kDa
Uniprot ID :	P42575

### Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IHC
Recommended dilution:	WB: 1:500-1:2000; IHC: 1:50-1:200
Image:	



Immunohistochemistry of paraffin-embedded rat kidney tissue slide using FNab01288(Caspase 2 Antibody) at dilution of 1:100



HeLa cells were subjected to SDS PAGE followed by western blot with FNab01288(CASP2 antibody) at dilution of 1:1000