

CBLL1 antibody

Product Information

Catalog No.:	FNab01320
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


Promotes ubiquitination of several tyrosine-phosphorylated Src substrates, including CDH1, CTTN and DOK1. Targets CDH1 for endocytosis and degradation(By similarity).

Immunogen information


Immunogen:	Cas-Br-M(murine) ecotropic retroviral transforming sequence-like 1
Synonyms:	E3 ubiquitin-protein ligase Hakai Casitas B-lineage lymphoma-transforming sequence-like protein 1 (c-Cbl-like protein 1) RING finger protein 188 RING-type E3 ubiquitin transferase Hakai CBLL1 HAKAI RNF188
Observed MW:	56 kDa
Uniprot ID :	Q75N03

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IHC, IF
Recommended dilution:	WB: 1:500-1:2000; IHC: 1:20-1:200; IF: 1:20-1:200
Image:	



Immunohistochemistry of paraffin-embedded human breast cancer using FNab01320(CBLL1 antibody) at dilution of 1:50



DU 145 cells were subjected to SDS PAGE followed by western blot with FNab01320(CBLL1 antibody) at dilution of 1:500