

## CBX4 antibody

### Product Information

Catalog No.:	FNab01331
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)


### Background

#### Immunogen information


Immunogen:	chromobox homolog 4 (Pc class homolog, Drosophila)
Synonyms:	E3 SUMO-protein ligase CBX4 Chromobox protein homolog 4 Polycomb 2 homolog (Pc2, hPc2) CBX4
Observed MW:	85 kDa
Uniprot ID :	O00257

#### Application

Reactivity:	Human, Mouse
Tested Application:	ELISA, IHC, IP, WB
Recommended dilution:	WB: 1:500 - 1:2000; IHC: 1:50 - 1:200; IP : 1:50 - 1:200
Image:	



Immunohistochemistry of paraffin-embedded human lung cancer using FNab01331(CBX4 antibody) at dilution of 1:100



Jurkat cells were subjected to SDS PAGE followed by western blot with FNab01331(CBX4 antibody) at dilution of 1:1000