

## CCDC115 antibody

### Product Information

Catalog No.:	FNab01345
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

### Background


May be involved in Golgi homeostasis.

### Immunogen information

Immunogen:	coiled-coil domain containing 115
Synonyms:	Coiled-coil domain-containing protein 115 CCDC115
Observed MW:	22-25 kDa
Uniprot ID :	Q96NT0

### Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB
Recommended dilution: WB:	1:500-1:2000
Image:	



mouse brain tissue were subjected to SDS PAGE followed by western blot with FNab01345(CCDC115 antibody) at dilution of 1:300