

NSRP1 antibody

Product Information

Catalog No.:	FNab01362
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


None

Immunogen information

Immunogen:	coiled-coil domain containing 55
Synonyms:	Nuclear speckle splicing regulatory protein 1 Coiled-coil domain-containing protein 55 Nuclear speckle-related protein 70 (NSrp70) NSRP1 CCDC55 NSRP70
Observed MW:	66 kDa70 kDa
Uniprot ID :	Q9H0G5

Application

Reactivity:	Human
Tested Application:	ELISA, WB
Recommended dilution: WB:	1:500-1:2000
Image:	



COLO 320 cells were subjected to SDS PAGE followed by western blot with FNab01362(CCDC55 Antibody) at dilution of 1:600