

COA3 antibody

Product Information

Catalog No.:	FNab01363
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


None

Immunogen information

Immunogen:	coiled-coil domain containing 56
Synonyms:	Cytochrome c oxidase assembly factor 3 homolog, mitochondrial Coiled-coil domain-containing protein 56 Mitochondrial translation regulation assembly intermediate of cytochrome c oxidase protein of 12 kDa COA3 CCDC56 MITRAC12
Observed MW:	12 kDa
Uniprot ID :	Q9Y2R0

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, IHC, IF
Recommended dilution:	IHC: 1:20-1:200; IF: 1:20-1:200
Image:	



Immunohistochemistry of paraffin-embedded human colon cancer slide using FNab01363 (CCDC56 Antibody) at dilution of 1:50