

## CCHCR1 antibody

### Product Information

Catalog No.:	FNab01375
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

### Background


May be a regulator of keratinocyte proliferation or differentiation.

### Immunogen information


Immunogen:	coiled-coil alpha-helical rod protein 1
Synonyms:	Coiled-coil alpha-helical rod protein 1 Alpha-helical coiled-coil rod protein Putative gene 8 protein (Pg8) CCHCR1 C6orf18 HCR
Observed MW:	85 kDa
Uniprot ID :	Q8TD31

### Application

Reactivity:	Human
Tested Application:	ELISA, WB, IHC, IF
Recommended dilution:	WB: 1:500-1:2000; IHC: 1:20-1:200; IF: 1:20-1:200
Image:	



Immunohistochemistry of paraffin-embedded human testis tissue slide using FNab01375(CCHCR1 Antibody) at dilution of 1:50



Transfected HEK-293 cells were subjected to SDS PAGE followed by western blot with FNab01375(CCHCR1 Antibody) at dilution of 1:1000