

## CCL27 antibody

### Product Information

Catalog No.:	FNab01382
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

### Background


Chemotactic factor that attracts skin-associated memory T-lymphocytes. May play a role in mediating homing of lymphocytes to cutaneous sites. Binds to CCR10.

### Immunogen information

Immunogen:	chemokine(C-C motif) ligand 27
Synonyms:	C-C motif chemokine 27 CC chemokine ILC Cutaneous T-cell-attracting chemokine (CTACK) ESkine IL-11 R-alpha-locus chemokine Skinkine Small-inducible cytokine A27 CCL27 ILC SCYA27
Observed MW:	None
Uniprot ID :	Q9Y4X3

### Application

Reactivity:	Human
Tested Application:	ELISA, IHC, IF
Recommended dilution:	IHC: 1:20-1:200; IF: 1:20-1:200
Image:	



Immunohistochemistry of paraffin-embedded human testis tissue FNab01382( CCL27 antibody) at dilution of 1:50