

CCL28 antibody

Product Information

Catalog No.:	FNab01383
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


Chemotactic activity for resting CD4, CD8 T-cells and eosinophils. Binds to CCR3 and CCR10 and induces calcium mobilization in a dose-dependent manner.

Immunogen information


Immunogen:	chemokine(C-C motif) ligand 28
Synonyms:	C-C motif chemokine 28 Mucosae-associated epithelial chemokine (MEC) Protein CCK1 Small-inducible cytokine A28 CCL28 SCYA28
Observed MW:	8 kDa-9 kDa
Uniprot ID :	Q9NRJ3

Application

Reactivity:	Human, Mouse
Tested Application:	ELISA, WB, IHC
Recommended dilution:	WB: 1:500-1:2000; IHC: 1:20-1:200
Image:	



Immunohistochemistry of paraffin-embedded human colon using FNab01383(CCL28 antibody) at dilution of 1:50



human spleen tissue were subjected to SDS PAGE followed by western blot with FNab01383(CCL28 antibody) at dilution of 1:300